

Solar Biologicals Inc. / Qualicum Scientific Ltd

Material Safety Data Sheet
Dipslide / Culture Media

EffectiveDate

Feb 1,2008

Section 1: Product and Company Information

Product No.	MC-500
Trade Name	Clinical Dipslide
Manufacturer	Solar Biologicals Inc. / Qualicum Scientific Ltd. 15 Capella Ct Unit 128 Ottawa, ON K2E 7X1 Tel 800-344-4652 Fax 888-246-7769 Email solar@solarbiologicals.com

Section2: Composition, Information on Ingredients

Description	CLED Agar w/Polymyxin B / EMB (Levine) Agar - Plastic Dipslide
-------------	---

Section 3: Hazards Identification

Most important hazards	Not classified as hazardous
------------------------	-----------------------------

Section 4: First Aid Measures

First Aid – Eyes	Wash out with plenty of water
First Aid – Skin	Wash skin with soap and water
First Aid – Ingestion	Wash out mouth and water. Have person drink 1-3 glasses of water to dilute stomach contents

Section 5: Fire Fighting Measures

Extinguishing media – Suitable	Not combustible
--------------------------------	-----------------

Section 6: Accidental Release Measures

Personal precautions	No specific measures necessary
Spill	Transfer into suitable containers for recovery or disposal. Finally flush area with plenty of water

Section 7: Handling and Storage

Handling	Avoid contact with eyes, skin and clothing
Storage	Storage temperature should be controlled between 8C – 22C

